



VIS & IR - Evidence Camera

for crime scene and laboratory



System Concept

SCENEview VisIR is a modular system which can be easily configured and optimized for each specific crime scene or evidence type.

The basic configuration of the system features an exceptionally bright camera that operates in both the visible and infrared spectra.

Attestors forensic lighting system, LIGHTcube, delivers modular, high-intensity illumination with exceptional performance across 12 wavelengths, ranging from ultraviolet to infrared.

Users can adjust the light intensity of each LIGHTcube, attach LIGHTcubes of the same color to further increase intensity, or combine different wavelengths to best suit the application.

Although SCENEview VisIR automatically adjusts the camera settings and selects the appropriate filter, users also have the option to manually control filter selection, focus, brightness, and contrast adjustments if preferred.

The basic setup includes:

- 12 MP VIS-IR colour camera with autofocus (8 cm to ∞)
- precision digital zoom from x1 to x10
- Large touchscreen display for intuitive operation
- Integrated control computer with photo and video recording capabilities – no additional PC required
- Interface for connecting LIGHTcubes
- Automated filter selection



very flexible and application-specific forensic illumination via LIGHTcubes

mounting point and control interface for LIGHTcubes (with swivel function for optimal lighting)

stable attachment for carrying strap

8" color display with touch function

shutter release for photo and video recording

VIS bandpass filter

USB & HDMI connectors

powerful battery with a typical runtime of >> 60 minutes

automatic filter selection system with manual override

3/8" screw connection for tripod

LIGHTcubes can also be operated independently of SCENEview via a separate handle. An optional camera adapter allows use with existing camera systems.



SCENEview VisIR

SCENEview VisIR is the latest advancement in our newly developed line of camera systems and accessories, designed specifically for crime scene and laboratory evidence photography.



This innovative system allows you to seamlessly search for, detect, and document various types of evidence, such as body fluids (including saliva, urine, and seminal fluid), blood, gunshot residues, fibres, hair, fingerprints, and more.

As the first system in our new forensic product line, SCENEview VisIR is equipped with a cutting-edge camera system featuring high light sensitivity, extending into the infrared (IR) spectrum. It offers fast autofocus and a user-friendly touchscreen interface, ensuring efficient and intuitive operation. The lightweight design focuses on essential forensic functionalities, removing unnecessary features while ensuring user comfort.

Upon activating a light source, SCENEview VisIR automatically optimises the camera settings and selects the appropriate filter for the specific application.

The system's adjustable light source angle enhances the working area and delivers superior imaging quality compared to systems with fixed multi-point LED illumination.

The image at the top right demonstrates how fingerprints on a plastic surface are visible using reflected light, though the surface texture of the plastic affects the image quality. Using a flatter angle of illumination, as shown in the bottom right image, significantly improves the clarity of the ridge detail.



SCENEview VisIR can also be used on a laboratory stand or on **PHOTOvent**.



Further components are in preparation.

Image Storage and Transfer

Integrated SSD hard drive with 32 GB
Image transfer via USB stick

Display

8" colour touch display with 1280 x 800 pixels

Image Formats

Maximum 4608 x 2592 pixels
formats JPG, TIF, BMP, PNG, DNG (RAW)

Dimensions (L x W x H)

202 x 253 x 130 mm
(without LIGHTcube)

Weight

Camera < 1 kg (without LIGHTcube)

Power Supply

Via battery (typical run time >> 60 minutes)
alternatively via mains power supply
or copy stand adaptor

Available Light Options

Daylight white (6500 K), neutral-white (5000 K), UV (365 nm), violet (410 nm), blue (447 nm), blue-green (490 nm), cyan (505 nm), green (530 nm), orange (590 nm), red (630 nm), IR I (850 nm), IR II (940 nm)



This is a product info brochure. Images might not be scale. Binding is solely the separately available te specification.
Attestor Forensics GmbH reserves the right to alter the or specification without prior notice.



Attestor Forensics GmbH
Bad Wurzach
Germany

attestor@attestor-forensics.com

www.attestor-forensics.com