

Dino-Lite AM5018MT

The Dino-Lite AM5018MT can be connected directly to a HDTV or a LCD monitor through a direct digital DVI interface. It offers sharp and vibrant HD 720p resolution to display even the smallest details. The AM5018MT is the ideal choice in situations, where real time images are vital, such as working on PCB's or other miniature objects. By tapping the Microtouch sensor on the microscope, the image can be frozen for inspection of important details and by pressing the Microtouch sensor for 2 seconds the LEDs can be switched on and off. The aluminium housing offers enhanced protection and durability and has the best look and feel.



Model	AM5018MT
Resolution	HD 720p
Magnification	10x ~ 70x, ~200x
Interface	DVI
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Housing	Full aluminium alloy body
Microtouch	√
Magnification lock	-
Video Frame Rate	Up to 60 fps
Filter/Diffuser	-
Other	Type C to Type A connector cable & Type C to DVI-D cable