

HD-M1

Brand: Dino-Lite
Product Code: HD-M1



Short Description

Aluminium holder for Dino-Lite

Description

High grade metal holder for all Dino-Lite models, compatible with 8mm and 10 mm pin connector stands from both Dino-Lite and third party vendors. This holder can be ordered separately but is also part of the MS-W1 package.

Specification

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	4.7cm (L) x 3.8cm (W) x 6.3cm (H)
Weight	150 g
Features	
ESD safe	No
Information	
Package contents	1 holder
Warranty information	2 years European warranty
Accessory range:	Holders and adapters
Price range	€0 - €50,00