

MS15X

Brand: Dino-Lite
Product Code: MS15X



Short Description

High precision X/Y table

Description

The MS15X is an advanced XY table that can be used independently or in combination with the RK-06A or RK-10A column stands.

The MS15X X/Y table base with detachable rotating table is designed for use with all sorts of applications that require fine examination in an X- and Y-axis movement. With the 360 degree rotating table, one can easily carry out all kinds of precision examinations. The MS15X can be used stand-alone or attached to the RK-06A and RK-10A stands (with the included C-Clamps)

Specification

Compatibility	
Compatible with:	All models; can be mounted on RK-06A/RK-10A
Housing	
Dimensions	12.7cm (W) x 15.2cm (L) x 7.6cm (H)
Weight	907 g
Features	
ESD safe	No
Information	

Package contents	Base, stage, clamps
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€150,00 - €250,00

Product Gallery

