

MS30A

Brand: Dino-Lite
Product Code: MS30A



Short Description

Tripod Adapter

Description

Adapter to use the Dino-Lite holder on third party tripods (standard thread).

Note: the tripod is not included.

Specification

Compatibility	
Compatible with:	HL2, HD-P1 and HD-M1 holder
Housing	
Dimensions	10mm (D) x 40mm (L)
Weight	24 g
Features	
ESD safe	No
Information	
Package contents	1 adapter
Warranty information	2 years European warranty
Accessory range:	Holders and adapters
Price range	€0 - €100,00

Product Gallery

