

MS53B

Brand: Dino-Lite
Product Code: MS53B

Short Description

Flexible arm stand

Description

Dino-Lite heavy duty jointed flex arm stand. The MS53B is a jointed articulating flex-arm type stand with holster for use with the Dino-Lite hand held digital microscopes.

The precision machine crafted components provide smooth adjustment and unlimited horizontal positioning in a well balanced ergonomic design. The MS53B can be easily clamped onto and secured to any workstation or table top. Manufactured from durable, lightweight recycled aluminum and able to be folded flat makes the versatile MS53B stand conveniently portable.

The available MS53BA2 option includes the MS53B flex-arm stand with an additional rigid floor stand and base assembly for use when a workstation table top is not convenient or available.

Specification

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	35.5cm (W) x 16.5cm (L) x 5cm (H)
Weight	1.8 kg
Features	
ESD safe	No
Information	
Package contents	Stand, HL2 holder
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€350,00 - €500,00

Product Gallery

