

N3C-D Diffuser Cap

Brand: Dino-Lite
Product Code: N3C-D



Short Description

Usage: Working with shiny objects (metal/plastics/PCB/jewelry etc.).

Description

The interchangeable front caps are designed to be used with Dino-Lite Edge and Edge^{PLUS} series models. The caps provide additional versatility to Dino-Lite digital microscope applications. This is a translucent cap that partially softens the LED light. When working with shiny objects (metal/plastics/PCB/jewelry, etc)

Compatibility

All Dino-Lite Edge/Edge^{PLUS} models

N3C caps overview

| | MODEL | VERSIO N | N3C-O | N3C-E | N3C-C | N3C-L | N3C- D/D2 | N3C-S |
|----|---|--|----------|----------------------|------------|----------|-----------------|------------------|
| | | | Open cap | Extended open cap | Closed cap | Long cap | Diffuser cap | Sidelight cap |
| AM | 4115, 4117, 4515, 4517, 4815, 491 | Standard models: ZT, ZTL, MZT, MZTL, | ? | ? | ? | ? | ? | ? |

| | | | | | | | | |
|----|---|---|-------------------------|---|---|---|---|---|
| | 7, 7115, 7515, 7915, 73115, 73515, 73915 | ZTW | | | | | | |
| AM | 4115, 4515, 4117, 4517, 7115 | Ultraviolet /Infrared models: FVW, FUW, FIT, FJT, FKT , FUT, FVT, | ? (black version) | x | x | x | x | x |
| AM | 4115, 4117, 4515, 4517, 4815 | Fluorescen t models: CFVW, DFRW, GFBW, G2FBW, GRFBY, YFGW, RFYW | ? (black version) | x | x | x | x | x |
| AM | 4515, 7515, 73515 | High magn ification models: ZT4, T5, T8, MT4A, MT8A | ? | x | ? | x | x | x |
| AM | 4115, 7115, 73115 | Extra Long Working distance models: ZT, ZTF | ? | x | x | x | x | x |
| AM | 5216, 5218 | VGA, HDMI models: ZT, ZTL | ? | ? | ? | ? | ? | ? |