

N3C-S Sidelight cap

Brand: Dino-Lite
Product Code: N3C-S



Short Description

Usage: Working with shiny objects (metal/plastics/PCB.jewelry etc.).

Description

The interchangeable front caps are designed to be used with Dino-Lite Edge and Edge^{PLUS} series models. The caps provide additional versatility to Dino-Lite digital microscope applications. Use the sidelight cap to light from the side for highlighting the texture and depth of targets at short working distances. When working with shiny objects (metal/plastics/PCB/jewelry,etc)

Compatibility

All Dino-Lite Edge/Edge^{PLUS} models

N3C caps overview

	MODEL	VERSIO N	N3C-O	N3C-E	N3C-C	N3C-L	N3C- D/D2	N3C-S
			Open cap	Extended open cap	Closed cap	Long cap	Diffuser cap	Sidelight cap
AM	4115, 4117, 4515, 4517, 4815, 491	Standard models: ZT, ZTL, MZT, MZTL,	?	?	?	?	?	?

	7, 7115, 7515, 7915, 73115, 73515, 73915	ZTW						
AM	4115, 4515, 4117, 4517, 7115	Ultraviolet /Infrared models: FVW, FUW, FIT, FJT, FKT , FUT, FVT,	? (black version)	x	x	x	x	x
AM	4115, 4117, 4515, 4517, 4815	Fluorescen t models: CFVW, DFRW, GFBW, G2FBW, GRFBY, YFGW, RFYW	? (black version)	x	x	x	x	x
AM	4515, 7515, 73515	High magn ification models: ZT4, T5, T8, MT4A, MT8A	?	x	?	x	x	x
AM	4115, 7115, 73115	Extra Long Working distance models: ZT, ZTF	?	x	x	x	x	x
AM	5216, 5218	VGA, HDMI models: ZT, ZTL	?	?	?	?	?	?