

RK-02

Brand: Dino-Lite

Product Code: RK-02



Short Description

Gooseneck stand

Description

The RK-02 is a durable gooseneck stand solution designed for flexible and versatile usage with a Dino-Lite digital microscope. The C-clamp of the RK-02 fits to virtually any surface of up to 7 cm (~3 inch) thick, and can be mounted horizontally or vertically. With this sturdy stand, you can easily position the Dino-Lite by bending the 60 cm long arm and adjusting the rotating head. The base of the RK-02 can also be rotated providing even more freedom and it can be fixed with the quick release fastener for use on uneven surfaces. The RK-02 gooseneck arm is the ideal choice to flexibly use any of the Dino-Lite models, except the high magnification models.

Specification

Compatibility	
Compatible with:	All Dino-Lite models, except High Magnification models
Housing	
Dimensions	Gooseneck 60cm (L) x 1.5cm (D)
Weight	1.5 kg
Features	
ESD safe	No

Information	
Package contents	Stand, HD-P1 holder
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€100,00 - €200,00

Product Gallery

