

RK-05



Brand: Dino-Lite
Product Code: RK-05

Short Description

Portable rack
Tilt and stay feature (30 degrees)
Fast assembly

Description

If you want to observe your object or specimen from a different perspective, the Dino-Lite RK-05 stand provides a multi-angle view with its 30 degrees inclinable pole. RK-05 is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-05F with fine focus adjustment](#)

Specification

Compatibility	
Compatible with:	All Dino-Lite models, except Long Working Distance models
Housing	

Dimensions	14cm (W) x 11cm (L) x 16.5cm (H)
Weight	611 g
Features	
ESD safe	No
Information	
Package contents	Stand, HD-P1 holder
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€150,00 - €250,00

Product Gallery

