

# RK-06A



**Brand:** Dino-Lite  
**Product Code:** RK-06A

## Short Description

Portable rack  
Fast assembly  
Fine-focus adjustment

## Description

The Dino-Lite RK-06A stand is a sturdy and stable mid-range stand solution that can be used with all Dino-Lite digital microscopes. It is constructed of stainless steel and lightweight aluminum and offers precise fine-focus adjustment as well as a quick vertical release function.

Its adjustable horizontal arm has a 15cm range with 360-degree rotation and with the RK-10-EX arm (included) its effective range can be extended to observe larger objects. Various optional accessories are available for the RK-06A, such as the RK-10-FX flexible arm or the RK-10-PX, XY arm. With these features, the RK-06A is the ideal mid-range stand solution for precise orientation and maximum flexibility.

[AVAILABLE EXTENSIONS FOR THE RK-06A STAND](#)

---

**RK-10-EX**



additional horizontal arm (included with RK-10A)

**RK-10-FX**



flexible arm extension

**RK-10-VX**



vertical arm extension

**RK-10-PX**



*(Temporarily unavailable)*

### **Specification**

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	22 cm (L) x 15 cm (W) x 27 cm (H)
Weight	1060 g
Features	
ESD safe	No

## Information

Package contents	Stand, HD-P1 holder, RK-10-EX horizontal arm
Warranty information	2 years European warranty
Accessory range:	Professional stands
Price range	€200,00 - €300,00

## Product Gallery

