

# RK-06F



**Brand:** Dino-Lite  
**Product Code:** RK-06F

## Short Description

Portable rack  
Fast assembly  
Fine-focus adjustment

## Description

The Dino-Lite RK-06F stand is a sturdy and stable mid-range stand solution, with fine focus adjustment that can be used with all Dino-Lite digital microscopes. It is constructed of stainless steel and lightweight aluminum and offers precise fine-focus adjustment as well as a quick vertical release function.

Its adjustable horizontal arm has a 15cm range with 360-degree rotation and with the RK-10-EX arm (excluded) its effective range can be extended to observe larger objects. Various optional accessories are available for the RK-06F, such as the RK-10-FX flexible arm or the RK-10-PX, XY arm. With these features, the RK-06F is the ideal mid-range stand solution for precise orientation and maximum flexibility.

## [AVAILABLE EXTENSIONS FOR THE RK-06A STAND](#)

---

### **RK-10-EX**



additional horizontal arm (included with RK-10A)

**RK-10-FX**



flexible arm extension

**RK-10-VX**



vertical arm extension

**RK-10-PX**



*(Temporarily unavailable)*

### **Specification**

|                      |                                   |
|----------------------|-----------------------------------|
| <b>Compatibility</b> |                                   |
| Compatible with:     | All Dino-Lite models              |
| <b>Housing</b>       |                                   |
| Dimensions           | 22 cm (L) x 15 cm (W) x 27 cm (H) |
| Weight               | 1060 g                            |
| <b>Features</b>      |                                   |
| ESD safe             | No                                |

| Information          |  |
|----------------------|--|
| Package contents     | Stand, HD-P1 holder, RK-10-EX horizontal arm |
| Warranty information | 2 years European warranty                    |
| Accessory range:     | Professional stands                          |
| Price range          | €250,00 - €400,00                            |

## Product Gallery

