

RK-06AE



Brand: Dino-Lite
Product Code: RK-06AE

Short Description

mid-range stand solution, ESD-safe
Precision stand with quick vertical release function
Various optional accessories are available

Description

Square metal base with ESD-safe coating with column and support providing several adaptable heights. The RK-06AE comes with an ESD (electrostatic discharge) safe base plate instead of the standard RK-06A base plate. The RK-06AE is a stable vertical desktop stand that was designed to be used with the complete line of Dino-Lite ESD-safe digital microscopes.

The Dino-Lite RK-06AE stand is a sturdy and stable mid-range stand solution that can be used with all Dino-Lite digital microscopes. It is constructed of stainless steel and lightweight aluminum and offers precise fine-focus adjustment as well as a quick vertical release function.

Its adjustable horizontal arm has a 15cm range with 360-degree rotation and with the RK-10-EX arm (included) its effective range can be extended to observe larger objects. Various optional accessories are available for the RK-06AE, such as the RK-10-FX flexible arm or the RK-10-PX, XY arm. With these features, the RK-06AE is the ideal mid-range stand solution for precise orientation and maximum flexibility.

AVAILABLE EXTENSIONS FOR THE RK-06A STAND

RK-10-EX



additional horizontal arm (included with RK-10A)

RK-10-FX



flexible arm extension

RK-10-VX



vertical arm extension

RK-10-PX



(Temporarily unavailable)

Specification

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	22 cm (L) x 15 cm (W) x 27 cm (H)
Weight	1060 g
Features	
ESD safe	Yes

Information	
Package contents	Stand, HD-P1 holder, RK-10-EX horizontal arm
Warranty information	2 years European warranty
Accessory range:	ESD stands
Price range	€200,00 - €300,00

Product Gallery

