

BL-ZW1

Brand: Dino-Lite
Product Code: BL-ZW1



Short Description

Polarized Backlight

Description

The BL-ZW1 backlight accessory opens up a whole new range of applications. The BL-ZW1 can be powered over USB or the DC adapter. The backlight has a built in, freely rotatable polarizer. The accessory is ideally suited to work with Dino-Lite models that have the polarization filter.

Specification

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	10.1cm (L) x 9.5cm (W) x 1.9cm (H)
Weight	330 g
Features	
ESD safe	No
Information	
Package contents	Backlight pack, USB cable, power adapter
Warranty information	2 years European warranty
Accessory range:	Light and control
Price range	€250,00 - €400,00

Product Gallery

