

SL-ZW1

Brand: Dino-Lite
Product Code: SL-ZW1



Short Description

Polarized Sidelight

Description

The SL-ZW1 is a polarized USB powered gooseneck light with the ability for the polarization to be adjusted by a simple rotation on the accessory itself. To take advantage of the polarizing light, the accessory is ideally suited to work with Dino-Lite models that have the polarization filter.

Specification

| | |
|----------------------|---------------------------|
| Compatibility | |
| Compatible with: | All Dino-Lite models |
| Housing | |
| Dimensions | 40cm (L) x 1.6cm (D) |
| Weight | 56 g |
| Features | |
| ESD safe | No |
| Information | |
| Package contents | 1 side light |
| Warranty information | 2 years European warranty |
| Accessory range: | Light and control |
| Price range | €0 - €100,00 |

Product Gallery

