

HDE-01



Brand: Dino-Lite
Product Code: HDE-01

Short Description

Vertical Dino-Lite handle

Description

The Dino-Lite HDE-01 handle is a lightweight aluminum handle which provides a more ergonomic and natural grip when using the Dino-Lite.

It is most commonly used as a grip for the Dino-Lite EarScope, but is also suitable for other applications where a vertical grip is more useful. The sturdy, lightweight aluminum makes this HDE-01 handle a good add-on for the Dino-Lite EarScope or any other Dino-Lite model.

Specification

Compatibility	
Compatible with:	All Dino-Lite models
Housing	
Dimensions	65mm base (W) x 166mm (H)
Weight	188 g
Features	
ESD safe	No
Information	
Package contents	1 holder, Quick start guide

Warranty information	2 years European warranty
Accessory range:	Holders and adapters
Price range	€100,00 - €200,00

Product Gallery



HDE-01 Net weight: 188 Grams (HDE-01 only)
Gross weight: 293 Grams (with packing box)