

DermaScope® Polarizer HR 200x

Brand: Dino-Lite Medical
Product Code: MEDL7DM



Short Description

5 Megapixel resolution
Robust metal housing
10-70x & 200x Magnification
Built-in polarizer
Medical Device Class 1

Description

The Dino-Lite DermaScope® Polarizer HR 200x (MEDL7DM) has a 5 megapixel camera and magnifications of 10 to 70x and around 200x. Also this model DermaScope® features a built-in and fully adjustable polarization filter to reduce the gloss effect of the skin. The MEDL7DM also stands out from the MEDL4DM by the sleek metal casing.

The Dino-Lite DermaScope® 200x models are special models at a much higher magnification, which is often used for specialized dermatological applications or, for example, the analysis of skin biopsies. Because of the high magnification of this model suspicious nevi can be analysed with more accuracy. In practice, this model is used for the analysis and interpretation of biopsies and other preparations, in short, more for laboratory research.



USB 2.0



Polarizer anti-reflection



Adjustable
~10-70x & 200x



5 megapixel
2592 x 1944



8 White LEDs
switchable



Metal
housing

Specification

Lighting	
Light/ LED type	White
Number of LEDs	8
LED on/off switchable:	Yes
Infrared filter	IR cut-filter >650 nm
Diffuser available	No
Polarizer	Yes, linear
Optics	
Magnification	10-70x, 200x
Lens type	Glass with anti-reflection coating
Sensor	
Sensor type	CMOS
Resolution	5 Megapixel (2592x1944)
Maximum frame rate	30 fps
Compatibility	
Interface	USB 2.0
Operating system	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Included software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Supported image formats (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Supported video formats (Windows)	WMV, FLV, SWF
Supported image formats (MacOS)	JPEG, PNG
Supported video formats (MacOS)	MOV (max 1.3MP)
Imaging standards	DirectShow, UVC
Housing	
Housing material	Aluminium alloy housing
Magnification lock	No
Dimensions	10.5 cm (L) x 3.2 cm (D)
Weight	140g
Cable length	1.8m

Features	
Special feature	No
Measurement	No
Calibration	No
Microtouch sensor	Yes
Information	
Package contents	Microscope, carry pouch, software CD, user manual
Warranty information	2 years European warranty
Regulatory approval	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745
Price range	€800,00 - €1000,00

Product Gallery

