

F25MM

Brand: Dino-Lite
Product Code: F25MM



Short Description

C-mount adapter for rigid endoscope

Description

C-mount adapter for rigid endoscopes with a focal length of 25mm. This C-mount adapter connects to rigid, straight through endoscopes with max. 32mm ocular diameter.

Specification

| | |
|----------------------|----------------------------|
| Compatibility | |
| Compatible with: | AM4023CT, AM4025X, AM7025X |
| Housing | |
| Dimensions | 48mm (H) x 52mm (W) |
| Weight | 120 g |
| Features | |
| ESD safe | No |
| Information | |
| Package contents | Adapter |
| Warranty information | 2 years European warranty |
| Accessory range: | Various |
| Price range | €250,00 - €350,00 |

Product Gallery

