

AF4515T-FUW

Designed for professionals in need of observing their target with different wavelengths, Dino-lite 1.3MP AF4515T-FUW digital microscope offers a 375nm UV and white light dual light source.



Overview

Specification

Model	AF4515T-FUW Dino-Lite Edge
Interface	USB 2.0
Product Resolution	1.3M pixels (1280x1024)
Magnification	20x~220x
Frame Rate	15fps in 1.3MP, MJPEG 30fps in VGA, MJPEG
Saving Format	Image: DinoCapture2.0: BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, RAS, PNM DinoXcope: PNG, JPEG Video: DinoCapture2.0: WMV, FLV, SWF DinoXcope: MOV
Lighting	4 white and 4 UV (375nm) LED lights switched by software
Microtouch	Yes
Operating System Supported	Windows XP/7/8/10/11 macOS 10.14 or later
Calibration Function	Yes
Measurement Function	Yes
Unit Dimension	11.4cm (H) x 3.3cm (D)
Unit Weight	111g