

# AM4115-FVT

The AM4115-FVT, with the exceptional image quality of Dino-Lite Edge featuring 395nm UV lighting with high sensitivity 1.3 Megapixels sensor., is designed to meet more professional needs. It also has better adaptability to different applications with the detachable and exchangeable front cap design.



## Overview

## Specification

Model	AM4115-FVT Dino-Lite Edge
Interface	USB 2.0
Product Resolution	1.3M pixels (SXGA)
Magnification Rate	20x~220x
Sensor	Color CMOS
Frame Rate	Up to 30fps
Save Formats	Image: DinoCapture2.0: BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, RAS, PNM DinoXcope: PNG, JPEG  Movie: DinoCapture2.0: WMV, FLV, SWF DinoXcope: MOV
Microtouch	Touch sensitive trigger on the microscope for taking pictures
Lighting	8 UV(395nm) LED lights switch on/off by software
Measurement Function	Yes
Calibration Function	Yes
Operating System Supported	Windows XP/7/8/10/11 macOS 10.14 or later
Unit Weight	100g
Unit Dimension	10.5cm (H) x 3.2cm (D)
Package Dimensions	16cm (L) x 16cm (W) x 6cm (H)