

# AM4115T-FVW

The AM4115T-FVW, with the exceptional image quality of Dino-Lite Edge, is designed to meet more professional needs. It also has better adaptability to different applications with the detachable and exchangeable front cap design.



## Overview

## Specification

|                            |  |
|----------------------------|--|
| Model                      | AM4115T-FVW Dino-Lite Edge   |
| Interface                  | USB 2.0  |
| Product Resolution         | 1.3M pixels (SXGA)   |
| Magnification Rate         | 20x~220x   |
| Sensor                     | Color CMOS   |
| Frame Rate                 | Up to 30fps  |
| Save Formats               | Image:<br>DinoCapture2.0: BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, RAS, PNM<br>DinoXcope: PNG, JPEG<br><br>Movie:<br>DinoCapture2.0: WMV, FLV, SWF<br>DinoXcope: MOV |
| Microtouch                 | Touch sensitive trigger on the microscope for taking pictures  |
| Lighting                   | 4 white and 4 UV(395nm) LED lights switched by software  |
| Measurement Function       | Yes  |
| Calibration Function       | Yes  |
| Operating System Supported | Windows XP/7/8/10/11<br>macOS 10.14 or later   |
| Unit Weight                | 100g   |
| Unit Dimension             | 10.5cm (H) x 3.2cm (D)   |
| Package Dimensions         | 16cm (L) x 16cm (W) x 6cm (H)  |