

# AM4115TL-FVW

The AM4115TL-FVW extends Dino-Lite's applications in long working distance to higher magnification range, featuring 395nm UV and white lights switchable by software with better image quality.



## Overview

## Specification

Model	AM4115TL-FVW Dino-Lite Edge
Interface	USB 2.0
Product Resolution	1.3M pixels (SXGA)
Magnification Rate	10x~140x
Sensor	Color CMOS
Frame Rate	Up to 30fps
Save Formats	Image: DinoCapture2.0: BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, RAS, PNM DinoXscope: PNG, JPEG  Movie: DinoCapture2.0: WMV, FLV, SWF DinoXscope: MOV
Microtouch	Touch sensitive trigger on the microscope for taking pictures
Lighting	4white and 4 UV(395nm) LED lights switched by software
Measurement Function	Yes
Calibration Function	Yes
Operating System Supported	Windows XP/7/8/10/11 macOS 10.14 or later
Unit Weight	100g
Unit Dimension	10.5cm (H) x 3.2cm (D)
Package Dimensions	16cm (L) x 16cm (W) x 6cm (H)