

Product data sheet

Dino-Lite AM4113ZT

By using a special polarization filter, the AM4113ZT is the right choice when working with high glare materials such as plastics and metals. It has a high resolution of 1,3 megapixel and an adjustable magnification of 10~70x as well as 200x. This model also provides the possibility to calibrate the microscope for measurements.



Model	AM4113ZT
Resolution	1,3 megapixel (1280x1024 pixels)
Magnification	10x~70x, 200x
Interface	USB 2.0, High speed USB transmission
Nr. of LED's	8
LED on/off switchable	√
LED color	White
Housing	Composite
Microtouch	√
Magnification lock	√
Output	Picture, video and time-lapsed video
Video Frame Rate	up to 30 frames-per-second, 15fps at 1.3M
Measurement	√ (line, radius, circle, 3-point circle, etc.)
Calibration	√
Filter/Diffuser	built-in polarizer
Compatibility	Windows XP, Vista, 7, 8 and recent MAC OS
Software (Included)	Windows: DinoCapture, Mac: DinoXcope
Other	The polarizer can be switched on/off or can be adjusted to offer full or half polarization effect

www.dino-lite.eu